

# Enterview® III - 5500

LED cordless examination headlight



The ENTERVIEW® III - 5500 is a rechargeable, cordless headlight which overcomes the need for connection to a fibre optic light source, separate battery pack or transformer, and hence increases freedom and mobility. The balanced design of the ENTERVIEW® III - 5500 headlight permits comfortable use during longer periods of time. The light emitting diode (LED) provides ample light for a regular outpatient session. The ENTERVIEW® III - 5500 is a versatile aid that is suitable for use in ENT, plastic surgery, gynaecology, maxillo facial surgery etc.

*Automatically switches on when taken off its battery charger. Switches off when placed back on its battery charger after use.*

## Benefits at a glance

- Cordless, rechargeable design
- Automatic on/off function
- Bright LED light
- Adjustable spot size
- Fully adjustable link between headband and lamp housing
- Adjustable headband with removable covering
- Multipurpose
- Approx battery capacity 3 hours
- Light weight

## Technical specifications

- Weight ca 325 grams
- Spot diameter (distance 25cm) adjustable 25-100mm
- Colour temperature 5500K
- Charging by specially designed battery charger
- Battery capacity sufficient for approx 3 hours use
- Luminous flux typical 235 lumen
- Collimating lens system
- Fresnel lens



Comfortable fit



Dimmer function



Adjustable spot size

### Product code: 1030220-WAND-ENT

Enterview® III – 5500 cordless headlight  
with wall-mounted battery charger

### Product code: 1030220-ENT

Enterview® III – 5500 cordless headlight  
with battery charger mounted on Entermed ENT treatment unit