

Morpheus Recline Cushion System

The Morpheus Recline Cushion System has been specifically developed for use with reclining chairs. It is designed to provide enhanced pressure distribution and relief to existing chair cushions, without compromising the patient's posture. The cushion can be easily expanded to include leg, back and arm sections providing high pressure reducing distribution and with the additional options for the Coccyx and Vertebral relief becomes a Very High Risk system offering offloading for the most vulnerable parts of the body.



Specifications

Product Code	Product Name	Dimensions	Weight	Risk Factor
①APH123CUSH	Main Cushion	47 x 43 x 5 cm	0.8 kg	High Risk
②APH123BACK	Back Cushion	47 x 70 x 3 cm	0.9 kg	
③APH123LEG	Leg Cushion	47 x 43 x 3 cm	0.5 kg	
④APH123ARMS	Arm Pads (pair)	20 x 43 x 2 cm	0.4 kg	
⑤APH123CUSHC R	Main Cushion with Coccyx Relief	47 x 43 x 5 cm	0.8 kg	Very High Risk
⑥APH123BACKV R	Back Cushion with Vertebral Relief	47 x 70 x 3 cm	0.9 kg	
APH123FULLWA	Complete Morpheus Recline Cushion System without Arms (①-③)			High Risk
APH123FULLA	Complete Morpheus Recline Cushion System with Arms (①-④)			
Maximum Weight Limit: 25 stone / 159 kg				
Risk Factor: High to Very High Risk				
Warranty: 1 year				
Safety Standards: BS EN 597-1 Cig, BS EN 597-2 Match, BS 7177 Ignition Source 5				